Get Energized! STEM Kit Clean Up

- 1. The multimeter needs to be turned off and the alligator clip cables should be unplugged and neatly coiled.
- 2. The digital balance and flashlight also need to be turned off.
- 3. The electrolyte solution is safe to be poured down the sink drain.
- 4. The electrolyte jar and ALL washers must be rinsed **thoroughly** and completely dried with a paper towel or soft cloth.

If the washers are not cleaned and dried properly before being put back into the bags, they will corrode, sticking together and be ruined.

- 5. All washers need to be placed into the appropriate bag (Solar A & Solar B copper washers go in the red bag all others go in the blue bag).
- Put the jar of salt inside of the dried, empty 100 mL beaker. Place the LED light, magnet and 4 circle papers into the empty jar with the weigh boat on top. <u>All items should be dry first!</u>



Get Energized! STEM Kit Packing Instructions

Step 1: Place the digital scale, multimeter, the red and blue cloth bags, and the 100 mL beaker with jar of salt rock in it into the case as shown.



Step 2: Put the flashlight on top of scale along with the jar with the magnet, LED & paper circles in it. The tweezers and pipette can be placed on top of the multimeter. The red and black cables should be gently wrapped up and put over the red and blue cloth bags.



Step 3: Lastly, place the Get Energized! booklet on top and carefully close the lid without pinching any of the items. <u>CAUTION:</u> *Do not force the case shut – it should close with ease!*

